

New Approach to 3-D Imaging of Biological Materials Provides Isotropic Resolution Down to a Few Nanometers

AL

Visualizing the three-dimensional (3-D) architecture of cells and tissues is essential for understanding the relationships between structure and function in biological systems. Biologists have primarily used light microscopy (LM) to investigate biological structures and examine living samples. Importantly, many biological samples are nearly transparent to visible light, permitting relatively easy 3-D analysis. However, the resolving power of LM is fundamentally limited by the wavelength of visible light: generally it cannot resolve details smaller than a few hundred nanometers.

Scanning electron microscopy (SEM) can resolve details down to the nanometer scale, sufficient for many of the structures of interest in current biological investigations. Unfortunately, 3-D analysis is challenging because biological materials are generally not transparent to electrons and tissues require heavy metal staining for contrast, necessary for SEM imaging.

Recent developments in a technique known as serial block face imaging (SBFI) allow automated acquisition and reconstruction of high-resolution, large tissue volumes.

Scanning electron microscopy and serial block face imaging

A scanning electron microscope scans a finely focused electron beam in a raster pattern over the surface of a bulk specimen and creates a virtual image by correlating changes in various signals with the position of the beam in the raster. The most useful signal for 3-D imaging is the

backscattered electron (BSE) signal, which reflects the average atomic number of the region illuminated by the beam. BSEs are beam electrons that are scattered back out of the sample surface by interactions with nuclei of sample atoms. They originate from a volume of interaction immediately below the beam, and the depth of that volume is determined by the energy of the beam. Biological electron microscopists have developed a number of staining techniques to enhance BSE contrast by selectively staining certain structures with heavy elements.

SEM image acquisition is a serial process, and requirements for resolution and signal-to-noise ratio (SNR) fundamentally determine the time required to acquire an image. Consequently, the strength of the signal and the efficiency are important considerations. The signal strength is determined mainly by the amount of beam current the sample can tolerate before damage occurs; biological samples, even when embedded, can be especially sensitive to beam damage. Therefore, the detection system needs to be optimized to collect low-loss BSE signal with high efficiency; this is best done using an in-lens detection geometry that generates images with strong contrast and high SNR.

Charging artifacts occur in SEM images when charge accumulates on non-conductive samples. Plastic-embedded biological samples can be very difficult to image. In routine imaging applications, conductive coatings may be applied to nonconductive samples to provide a path to ground that prevents charge accumulation. Coatings can be used in SBFI after a fresh surface is exposed with an in situ microtome; however, this increases




Figure 1 – The serial block face imaging workflow involves chemically fixed, heavy metal stained and resin-embedded biological sample block attached on an Al stub with conductive glue. Several regions of interest can be selected on the block face and imaged at high resolution with either individual fields of view or tile sets stitched together to form a large composite image. A layer of material is removed using a diamond knife, and the multi-energy acquisition (MEA) process is repeated, resulting in a stack of images that can be reconstructed by a computer into a high-resolution representation of the sectioned volume. Segmentation of the data can be used to isolate structures of interest.

the experimental time substantially. Low vacuum (LoVac) operation provides a means to eliminate charging that is compatible with SBF. In LoVac mode, a gas, in this case water vapor, is injected into the sample chamber and maintained at a pressure low enough to permit the beam to reach the sample without excessive scattering. Some molecules of the gas are ionized by beam electrons or by electrons originating from the sample. These ionized gas molecules are then available to neutralize any charge that accumulates on the sample. LoVac operation is essential to ensure high-quality imaging of difficult samples.

SBFI generates a sequence of images of the sample surface as the sample volume is serially sectioned, each section removing a thin layer of material from the block face (see Figure 1). The resulting stack of images can be computationally reconstructed to visualize the original 3-D structure. In an SEM, the lateral resolution of the reconstruction can be as good as a few nanometers; however, the depth resolution is limited by the thickness of the slices removed from the surface. This can be quite thin, since the structural integrity of the removed




Figure 2 – Multi-energy acquisition captures several images at different beam energies (landing energies). Each image captures information from a different depth below the block face surface, within the next layer of material to be removed. A proprietary software algorithm deconvolutes the data to provide a subset of optical sections for each physical section. The final reconstructed model exhibits isotropic resolution down to a few nanometers.

If you have your own static problem, NRD has the SOLUTION!

Did you know that even a small static charge can affect accuracy when weighing fine powders and particulates?

Static on weigh boats, filter paper - even clothing can cause inaccurate measurements. Eliminate the problem of static charge with NRD's complete line of anti-static devices designed specifically for precision weighing applications.

In chamber static elimination for fine powders and small particulates

Wash top static elimination for irregularly shaped parts

20500 Ionizing Cartridge with 5000 Counts

10 Min. Use

Contact the experts at NRD to find the solution that's right for your lab.

Digital Water Baths at analog prices.

Large, digital display shows set and actual temperature simultaneously.

Clear, hinged cover tilts out of the way and lets condensation drain back into the bath, while keeping samples in view.

Handles make for easy positioning and relocation of the bath.

Built-in timer eliminates need for an extra device.

Reach set temperature in a fraction of the time.*

Three user-defined temperature presets can be used for frequently used set temperatures.

USES 1/3* LESS ELECTRICITY

See us at Pittcon - Booth # 4127

PolyScience®
Temperature Control Solutions
www.polyscience.com
1.800.229.7569

Circulating Baths
General Purpose Water Baths
Chillers
Specialty Products
*Compared to previous models.




Figure 5 – Volume reconstruction of a mouse brain volume. a) A representative SEM image. b) The segmented 3-D reconstruction has isotropic resolution (x , y and z) of 10 nm. The physical sections were 50 nm thick. A total of 1040 sections, physical and optical, were used to reconstruct the $15.00 \mu\text{m} \times 12.9 \mu\text{m} \times 10.4 \mu\text{m}$ volume.

Mitochondrial network

Figure 5 shows the reconstruction of a complex mitochondrial network within neurons of mouse brain. Mitochondria are highly dynamic organelles that have central importance for ATP production and undergo fusion and fission events continuously, leading to a

diverse array of mitochondrial morphologies that range from fragmented states to continuous networks. Three-dimensional volume imaging allows reconstruction, visualization and analysis of these highly complex organelles within the context of the biological tissue under investigation.

Conclusion

A new approach to 3-D volume imaging of biological samples combines mechanical sectioning with optical sectioning using multi-energy deconvolution to provide isotropic resolution. Optimized detectors and low-vacuum operation ensure high-quality imaging. Software provides automation of processes ranging from low-level setup and alignment to walkaway acquisition of a complete image series. The approach includes a complete workflow from initial setup to final results, including software solutions for large area/volume analysis, LM image overlay, reconstruction, visualization and segmentation.

Emine Korkmaz is product marketing manager, Cellular Biology, FEI, Achtseweg Noord 5, 5651 GG Eindhoven, The Netherlands; tel.: +31 40 235 6000; e-mail: Emine.Korkmaz@fei.com; www.fei.com


The World's Most Reliable Equipment



The 914 semi-automatic changeover manifold available in brass and stainless steel is a workhorse in many processes and laboratories.



SGD protocol stations make changing cylinders easier while protecting your regulator from damage.



The AUTO-LOGIC II® touch screen electronic automatic changeover manifold provides high technology control at a low technology cost. Available in brass and stainless steel.



SGD point of use panels are a convenient organized way to deliver gases to your instrumentation or processes.



SGD provides a comprehensive line of pressure regulators for virtually any application.



Whisper Valve®
Takes the **POP!** out of cryogenic safeties and pays for itself in gas savings.

- Relieves cryogenic container pressure at less than 40 dB. Normal safeties pop at over 100 dB (the sound of a shotgun blast).
- Reduces container losses to 2 SCFH (maintains liquid in a cryogenic container for up to 50 days).
- Easily installed on any cryogenic container.
- Available in a range of pressure settings.



The improved design 8010 and 8000 cartridge gas purifiers remove a variety of impurities from gas streams at a reasonable cost.



23-27 Bland Street
Emerson, NJ 07630-0347
Phone: 201.261.1908
Fax: 201.265.2162
sgd@sgd.com
www.sgd.com